

Gaussian Beam Optics

Determining the intensity profile of a refracted laser beam

Louis Claeys, Ruben Van der Borght

Kevin Anthony Kaw
KU Leuven

May 23, 2022

Theory




Figure 1: The focal point of a thin lens¹.




Figure 2: A beam and its embedded Gaussian beam².

¹M. Run, Thin Lens Formula, Wikipedia (2019)

²CVI Melles Griot, Gaussian Beam Optics, Technical Guide 2009, 2, 1 (2009)

Theory




Figure 3: The radial Gaussian intensity profile $I_z(w)$ of TEM_{00}^3 .




Figure 4: A beam and its embedded Gaussian beam³.

³CVI Melles Griot, Gaussian Beam Optics, Technical Guide 2009, 2, 1 (2009)

Theory

- $I_z(w)$ is Gaussian for all z for TEM₀₀.
 - $w(z) = w_0 \sqrt{1 + \left(\frac{z\lambda M^2}{\pi w_0^2}\right)^2}$
-
- M^2 quality factor
 - Gaussian Beam: $M^2 = 1$
 - Real laser: $1.1 < M^2 < 1.3$
 - Study of non-linear absorption of clusters⁴

⁴Su Z. et al, Nanoscale, 14, 3618. (2022)

Aim

- Check if our laser beam is Gaussian.
- Determine M^2 for two wavelengths.

Setup




Figure 5: Schematic representation of the experimental setup⁵.




Figure 6: The experimental setup including the camera on the motorized linear stage.

⁵Su Z. et al, *Nanoscale*, 14, 3618. (2022)

Laser

- Quanta-Ray pulsed Nd:YAG (532 nm) pumped Sirah Cobra Stretch dye laser system
- $\lambda_1 = 600 \text{ nm}$, $\lambda_2 = 570 \text{ nm}$




Figure 7: The beam profile.



Figure 8: The used Nd:YAG laser (left) and the dye laser (right).

Camera




Figure 9: Thorlabs BC106N Beam Profiler.





Figure 10: Thorlabs Beam 7.7 (Beam Analyzing Software).

Measuring the profile: fitting Gaussians



$$I = Ae^{\frac{(z-\mu)^2}{2\sigma^2}}$$

$$w(z) = 2\sigma(z)$$

Figure 11: Gaussian profile data of the beam for $z = 17 \text{ cm}$ in the x-direction, along with fit.

Measuring the profile: error analysis $I(x), I(y)$




Figure 12: Standard deviation in repeatedly measured data for $z = 9 \text{ mm}$ ($n = 5$).

Measuring the profile: error analysis $w(z)$




Figure 13: Difference between w measured in the x and y direction (570 nm).




Figure 14: Difference between w measured in the x and y direction (600 nm).

Measuring the profile: fitting $w(z)$ ($\lambda = 600$ nm)




Figure 15: Measured laser profile along with fit for $\lambda = 600$ nm. $\chi^2_{\text{red}} = 1.18$, $p = 0.265$.

$$w(z) = w_0 \sqrt{1 + \left(\frac{z - z_0}{z_R} \right)^2}$$

$$w_0 = (13.0 \pm 0.3) \mu\text{m}$$

$$z_0 = (11.9 \pm 0.1) \text{ mm}$$

$$z_R = (0.81 \pm 0.02) \text{ mm}$$

$$M^2 = \frac{\pi w_0^2}{z_R \lambda} = 1.09 \pm 0.06$$

Measuring the profile: fitting $w(z)$ ($\lambda = 570 \text{ nm}$)




Figure 16: Measured laser profile along with fit for $\lambda = 570 \text{ nm}$. $\chi^2_{\text{red}} = 1.35$, $p = 0.149$.

$$w(z) = w_0 \sqrt{1 + \left(\frac{z - z_0}{z_R} \right)^2}$$

$$w_0 = (12.7 \pm 0.3) \mu\text{m}$$

$$z_0 = (10.0 \pm 0.1) \text{ mm}$$

$$z_R = (0.82 \pm 0.02) \text{ mm}$$

$$M^2 = \frac{\pi w_0^2}{z_R \lambda} = 1.08 \pm 0.06$$

Discussion: Chromatic Aberration




Figure 17: Data and fits for both wavelengths together: notice the shape is nearly identical, but the focal point has shifted due to chromatic aberration.

Conclusion

1. Gaussian beam approximation holds
2. $w(z)$ behaviour corresponds with theory ($p = 0.265, p = 0.149$)
3. $M^2(600 \text{ nm}) = 1.09 \pm 0.06, M^2(570 \text{ nm}) = 1.08 \pm 0.06$
4. Compatible with each other and suited for high precision experiments
5. Main source of errors: angular distribution

References

[Run, 2019] Run Mike (2019)

Thin Lens Formula

Wikipedia, https://commons.wikimedia.org/wiki/File:Thin-lens-formula-part1_DE.svg.

[CVI, 2009] CVI Melles Griot (2009)

Gaussian Beam Optics

Technical Guide 2009, 2, 1.

[Su, 2022] Su Z. et al. (2022)

The wavelength-dependent non-linear absorption and refraction of Au_{25} and Au_{38} monolayerprotected clusters

Nanoscale 2022, 14, 3618.

Methods

1. Align the camera.
2. Adjust z to the desired value.
3. Uncover the laser and calibrate the camera.
4. Aim the crosshair on the peak value.
5. Collect the data.
6. Repeat for several values of z .





Figure 18: The camera alignment and software.